



## **BrainstormTechnion Hack Nights 2021**

### **Intro to EEG & OpenBCI workshop**

**Or Rabani**

#### **Getting started with OpenBCI hardware**

1. Download openBCI GUI version 4.2 from their website.  
This version has the Focus Widget.
2. Connect your electrodes (1&2) to the board using this manual:  
If your board is white:  
<https://docs.openbci.com/GettingStarted/Boards/CytonGS/>  
If your board is blue:  
<https://docs.openbci.com/GettingStarted/Boards/GanglionGS/>
3. Use Gel and ribbon to connect the electrodes to your forehead  
follow this guide and the instruction I gave at the workshop  
<https://docs.openbci.com/GettingStarted/Biosensing-Setups/EEGSetup/>
4. In the GUI, open 4 window layout (top right), filter at 50HZ and BP1-50 (top left)
5. Open: raw signal, FFT, Networking and impedance widgets.
6. Disable channels 3 & 4 by clicking on the colored circle with the number
7. check for impedance using the impedance weight
8. replace Impedance widget with focus widget.  
<https://docs.openbci.com/Software/OpenBCISoftware/GUIWidgets/>
9. See if you can identify artifacts like eye blinks, and muscle activity
10. Focus and have fun :)

## BrainstormTechnion Hack Nights 2021

### Intro to EEG & OpenBCI workshop

Or Rabani

#### Getting started with OpenBCI GUI

1. Download openBCI GUI version 4.2 from their website.  
This version has the Focus Widget.
2. Download and install the latest Arduino IDE  
<https://www.arduino.cc/en/software>
3. Download the \*.ino file from step 4 in this tutorial  
<https://www.instructables.com/Send-Focus-Data-From-OpenBCI-GUI-to-Arduino/>
4. Connect the Arduino and upload the file ("sketch") to the board using Arduino IDE and this manual:  
<https://www.arduino.cc/en/Tutorial/getting-started-with-ide-v2/ide-v2-uploading-a-sketch>
5. Connect an LED to pin 13 (long leg) and the GND pin next to it (short leg)
6. Load Or's brain from here:  
[https://drive.google.com/file/d/1Om3J8KliprQp6pGm9EV0b5Uxa\\_1bhD8M/view?usp=sharing](https://drive.google.com/file/d/1Om3J8KliprQp6pGm9EV0b5Uxa_1bhD8M/view?usp=sharing)  
In the GUI: choose load from file, 4 channels.  
Use this manual to get to know the GUI  
<https://docs.openbci.com/Software/OpenBCISoftware/GUIDocs/>
7. In the GUI, open 4 window layout (top right), filter at 50HZ and BP1-50 (top left)
8. Disable channels 3 & 4 by clicking on the colored circle with the number
9. Open: raw signal, FFT, Networking and focus widgets.  
More about the focus widget:  
<https://docs.openbci.com/Software/OpenBCISoftware/GUIWidgets/>
10. connect between Arduino and OpenBCI using steps 6-7 in this manual  
<https://www.instructables.com/Send-Focus-Data-From-OpenBCI-GUI-to-Arduino/>  
See if you can identify when I am blinking, moving or focusing :)